

Multi-cycle high-throughput growth media optimization using batch Bayesian optimization

Frédéric M. Lapierre^{1,2}, Pasquale Mattaliano², Dennis Raith^{2,3}, Mariela Castillo-Cota^{2,3}, Jonas Bermeitinger², Neda Sabet², and Robert Huber¹

Abstract

Optimizing microbial media is vital for industrial yields, yet traditional **Design of Experiments (DoE)** often **misses complex interactions**. This study evaluates **Batch Bayesian Optimization (BBO)** using *Sporosarcina pasteurii*.

In-silico, BBO outperformed DoE by **14%**. In automated **wet-lab trials**, BBO achieved a **28%** increase in biomass titer via **superior exploration of media components**. These results scaled successfully to a **2L bioreactor**, where BBO yielded **19%** more biomass than DoE. Overall, BBO provides a **more efficient, targeted** framework for **enhancing bioprocess yields** and **reducing costs**.

Introduction: Navigating Optimization Methods

Industrial fermentation success hinges on cost-effective, strain-specific culture media. Traditional **one-factor-at-a-time (OFAT)** methods are inefficient and miss complex interactions, while **Design of Experiments (DoE)** often struggles with the nonlinear nature of microbial growth. Bayesian Optimization (BO) offers a superior alternative by using Gaussian process surrogate models to balance exploration and exploitation. Specifically, **Batch Bayesian Optimization (BBO)** enables the simultaneous identification of multiple sampling points, drastically reducing the time required for optimization cycles in slow-growing biological systems.




Figure 1. Different optimization approaches: **a.** One-factor-at-a-time, which focuses on changing variables one by one, neglecting the interaction of different factors. **b.** Design of Experiments, where multiple factors are optimized together, and the interaction of factors is considered partially, and **c.** Automated Batch Bayesian Optimization, where multiple factors are optimized simultaneously, and complex interactions of factors are taken into account.

Methodology: The BBO Framework & Automated Hardware

The methodology employed in this study combines advanced **probabilistic modeling** with **high-throughput automation** to optimize the growth of *Sporosarcina pasteurii*. The BBO process begins by constructing a **Gaussian process (GP) model** from an initial data set, which provides both **function estimates** and **uncertainty measurements** critical for navigating the search space. To manage the computational complexity of selecting multiple points for a batch, the system adopts a heuristic strategy using local penalization to define exclusion zones around selected points, ensuring diverse exploration without redundant evaluations.

This study optimized the growth of *Sporosarcina pasteurii* by integrating **probabilistic modeling** with high-throughput automation. The BBO framework utilizes Gaussian processes to map the search space and local penalization heuristics to ensure **diverse, non-redundant** batch selection. This algorithmic framework is integrated with a **fully automated microbioreactor setup** consisting of a BioLector I system for **parallel cultivation** and **real-time monitoring** of biomass via backscatter. A **liquid handling system** (OT-2) is utilized to dispense media components directly into the microplates, while a small-scale turbidostat system ensures the continuous provision of a consistent inoculum for multiple automated cycles.





Figure 2. a. Automated cultivation device including BioLector I system, liquid handling system (OT-2), and turbidostat system. Right panel (b): Growth data for the microbioreactor cultivation of bacterium *Sporosarcina pasteurii*.

Results: Superior Exploration and Scalability

BBO outperformed DoE by maintaining a broader exploration of the design space. While DoE dismissed iron early as insignificant, BBO's iterative approach identified it as a critical limiting factor at high biomass concentrations. Consequently, the BBO-optimized medium achieved a **14% increase in-silico simulation**, **28% increase** in maximum backscatter in **microscale trials**. In a **2L laboratory bioreactor**, BBO demonstrated excellent scalability, yielding an OD600 of 67.4, a **19% improvement** over DoE and the highest batch cultivation density for *S. pasteurii* reported in literature to date.




DoE vs. BBO variable exploration and optimization results. As shown above, BBO explores the factors to a greater extent than DoE, which explores only a few values for each component. The graph also shows BBO's superior performance over the same number of optimization runs.

Outlook

Future investigations should **expand** the application of BBO to optimize actual **target product concentrations** and **fed-batch process parameters** like feed composition and rates. Additionally, the effectiveness of BBO should be further assessed against other advanced optimization algorithms and applied to a **broader range of objectives**, such as **optimizing product yield per media cost**.

Another further investigation following this study is **integration of BBO with high-throughput automation** which marks a significant transition toward fully autonomous, **set-and-forget intelligent systems** in bioprocess development. This **closed-loop learning system** can **automate the entire experimental lifecycle**; from **high-precision liquid handling** and **media formulation** to **real-time imaging readout**, **AI-based analysis**, and **optimization**.



High-throughput self-optimizing system including plate handler, microscope, pipetting robot, and incubator at the hardware level, equipped with experiment control, automated data acquisition, AI analysis, and BBO optimization at the software level.



1. Department of Engineering and Management, HM Hochschule München University of Applied Sciences, Munich, Germany
2. Labmaite GmbH, Freiburg im Breisgau, Germany
3. Department of Computer Science, University of Freiburg, Freiburg im Breisgau, Germany



info@labmaite.com

labmaite.com

Türkheimerstraße 4a, 79110 Freiburg, Germany

